

RNA FISH Protocol for Adherent Cells

I . **Product Description**

RNA FISH Probes consist of 25-48 fluorescently labeled oligonucleotides specifically designed to bind to RNA or DNA targets. These probes enable the detection of single RNA molecules as diffraction-limited spots using conventional or confocal fluorescence microscopy.

II . **Content of the Kit**

Name	Cat#	Storage
FISH probe		-80°C
10x Hybridization buffer	TDP-HB-25	-80°C
5x Washing buffer	TDP-WB-25	4°C
20x SSC	TDP-SSC-25	Room Temperature
Positive control		-80°C

III . **Materials without in the kit**

Name	Vendor	Cat.#
16% paraformaldehyde solution	Electron Microscopy Sciences	15710
10X Phosphate Buffered Saline (PBS)	Fisher Scientific	SH3002803
Ethanol for molecular biology	Fisher Scientific	04-355-450
Nuclease-free water	Thermo Fisher Scientific	10-977-015
Deionized formamide	Millipore	4650-500ML
4',6-diamidino-2-phenylindole (DAPI)	Thermo Fisher Scientific	D3571
20x SSC buffer	Invitrogen	AM9770
Mounting Medium	Thermo Fisher Scientific	P10144
12 mm round coverslips	California Bio-supply	633029

IV . Other Reagents and Consumables:

1. 12- or 24-well culture plates
2. RNase-free consumables (e.g., pipette tips)
3. Humidified chamber (or equivalent): 150 mm tissue culture plate or plastic box with a flat bottom, lined with water-saturated paper towels and a single layer of Parafilm® placed on top
4. Super-frost Plus Microscope slides
5. 37°C laboratory oven

V . Buffer Preparations

1. Fixation Buffer:

The final composition should be 4% (vol./vol.) formaldehyde in 1X PBS.

For a final volume of 10 mL, mix:

- 2.5 mL 16% formaldehyde solution
- 7.5 mL 1X PBS

2. Hybridization Buffer:

The final composition should be 15% (vol./vol.) formamide in Hybridization Buffer.

Prepare fresh Hybridization Buffer for each experiment. To account for solution viscosity, it is recommended to prepare an excess of 5% in volume.

For a final volume of 1 mL, mix:

- 850 µL RNA FISH Hybridization Buffer (9 mL)
 - 150 µL deionized formamide
- Mix thoroughly by vortexing and pipetting up and down.

3. 1x Wash Buffer:

Prepare the 1x Wash Buffer freshly for each experiment. For a final volume of 10 mL, mix:

- 2 mL 5x TriDixBio RNA FISH Wash Buffer
- 6.5 mL nuclease-free water

- 1.5 mL deionized formamide
Mix thoroughly by gentle vortexing.

4. DAPI Nuclear Stain Solution:

- 4',6-diamidino-2-phenylindole (DAPI) dissolved in 1x Wash Buffer at 5 ug/mL.
This solution is used for DAPI staining for use after hybridization.

5. Mounting Media:

- Prolong mounting media: Use 10 µL for each coverslip.

VI . Protocol for Adherent Cell Fixation

1. Fixation of Adherent Cell Lines

- 1) Grow cells on 12 mm round coverslips in a 12-well or 24-well cell culture plate.
- 2) Aspirate the growth medium and wash cells once with 0.5 mL (for 24-well) or 1 mL (for 12-well) of 1x PBS.
- 3) Add 0.5 mL (for 24-well) or 1 mL (for 12-well) of fixation buffer (4% PFA in 1x PBS).
- 4) Incubate at room temperature for 10 minutes.
- 5) Wash twice with 0.5 mL (for 24-well) or 1 mL (for 12-well) of 1X PBS.
- 6) To permeabilize, immerse cells in 1 mL of 70% (vol./vol.) ethanol for at least 2 hours at 4°C.
- 7) Cells can be stored at 4°C in 70% ethanol for up to one week before hybridization.

2. Hybridization, Washing, and Sealing Coverslips in Adherent Cells

- 1) If frozen before use, warm the reconstituted probe solution to room temperature. Mix well by vortexing and briefly centrifuge.
- 2) To prepare the Hybridization Buffer containing the probe, add 1 µL of probe stock solution to 50 µL of Hybridization Buffer with 15% formamide. Vortex and centrifuge. This creates a working probe solution of 1 ng/µL, enough for one coverslip, which will be used in steps 6 and 7.
- 3) Aspirate the 70% ethanol from the coverslip containing adherent cells within the 12-well or 24-well plate.

- 4) Add 1 mL of 2x SSC buffer with 15% formamide and incubate at room temperature for 5 minutes.
- 5) Prepare the humidified chamber by placing water-saturated paper on the bottom, and insert a flat rack wrapped with Parafilm into the chamber.
- 6) Within the humidified chamber, dispense 50 μ L of the Hybridization Buffer containing the probe onto the Parafilm on top of the rack.
- 7) Gently transfer and flip over the coverslip, with cells side down, onto the 50 μ L drop of Hybridization Buffer containing the probe.
- 8) Cover the humidified chamber with the lid and wrap it with foil to protect from light.
- 9) Incubate the chamber in the dark at 37°C for at least 4 hours (incubation can be extended up to 16 hours).
- 10) Gently transfer the coverslip, with cells side up, to a fresh 12-well plate containing 1 mL of 1x Wash Buffer.
- 11) Incubate in the dark at 37°C for 30 minutes.
- 12) Aspirate the 1x Wash Buffer, then add DAPI (final concentration of 5 ng/mL) within 1 mL fresh 1x Wash Buffer for nuclear staining and second washing.
- 13) Incubate in the dark at 37°C for another 30 minutes.
- 14) Aspirate the DAPI staining Wash Buffer, then add 1 mL of 2x SSC buffer. Incubate at room temperature for 5 minutes.
- 15) Add a small drop (~10 μ L) of Prolong mounting media to a glass slide. Use Kim-wipes paper to soak up 2x SSC buffer from the edge of coverslip. Mount the coverslip onto the slide with cells side down.
- 16) Gently wick away any excess mounting media from the perimeter of the coverslip (optional).
- 17) Place the sealed slide into a drawer at room temperature to avoid light for at least 3 hours or overnight.
- 18) The slides with sealed coverslips are ready for imaging.

3. Imaging

1) Microscope Requirements:

- Wide-field fluorescence microscope (e.g., Leica TIRF & D-STORM/THUNDER) or Confocal (e.g., Confocal SP8 Light-Sheet)
- High numerical aperture (>1.3) and 60-100x oil-immersion objective
- Strong light source (e.g., mercury or metal-halide lamp; newer LED-based light sources may also suffice)
- Appropriate filter sets for fluorophores

- Standard cooled CCD or sCMOS camera (ideally optimized for low-light-level imaging with a pixel size of 13 μm or less)